3135

Validation of a bioinformatic model for classifying non-tumor variants in a cell-free DNA liquid biopsy assay

Jennifer Yen¹, Pegah Safabakhsh¹, Jeffrey L. Werbin¹, Denis Tolkunov¹, Yu Fu¹, Katie Quinn¹, Robert Foley¹, Minh Tran¹, Xianxian Liu¹, Carin Espenschied¹, Scott Higdon¹, Brett Kennedy¹, Shile Zhang¹, Han-Yu Chuang¹ ¹ Guardant Health, Palo Alto, CA

Introduction

Liquid biopsy is a powerful, non-invasive tool for profiling tumors and identifying clinically relevant variants.

The presence of clonal hematopoiesis (CH) variants, and biological noise, due to aging and therapy has potential to confound biomarker interpretation.

Currently, comprehensive methods to filter out non-tumor variants require genotyping the white blood cell (WBC) fraction of the paired plasma sample, which is a costly, complicated workflow.

A plasma-only, bioinformatics solution to identify non-tumor variants is needed for accurate biomarker assessments in the cell-free DNA (cfDNA).

Methods

Variant calls were obtained from >250,000 plasma samples comprising healthy donor, early and late-stage cancer patients sequenced on the Guardant360TM GuardantREVEALTM. GuardantOMNI[™] and GuardantInfinity[™] liquid biopsy panels as well as public tissue datasets.

The model was trained on paired plasma and WBC datasets and optimized with 10-fold cross-validation to produce a non-tumor and tumor variant classifier.

To validate these calls, an independent cohort of 72 paired plasma and WBC advanced cancer samples genotyped on the were GuardantInfinity[™] assay. A cohort of 76 healthy genotyped on the samples, donor GuardantOMNI assay was also assessed.

Input Data



Figure 1. Model design. Features were engineered from Guardant internal and external public datasets and trained using 10-fold cross validation using multiple models. Only results from the Logistic Regression model are shown. Model validation was performed on an independent cohort of paired plasma and WBC late-stage samples sequenced on the GuardantInfinityTM panel, and healthy donors sequenced on the GuardantOMNI[™] panel.

Model performance: High ROC AUC and accuracy for predicted calls

A. Validation: Late-stage cohort of 72 paired plasma and WBC samples



B. Validation: Healthy donor cohort of 76 paired plasma and WBC samples

All somatic SNV/Indels in healthy donors are considered non-tumor



Figure 2. Model performance. Predictions for tumor and non-tumor status were compared to WBC confirmation in A) 713 somatic SNV/Indels from 72 paired plasma and WBC GuardantInfinity[™] samples and B) 243 somatic SNV/Indels from 76 paired plasma and healthy donors on GuardantOMNI[™]. Lower confirmation rate in WBC sequencing observed for low VAF variants (<0.6%) likely attributed to the limit of detection in for WBC variant calling and/or possible non-WBC lineage origin.



		Non-Tu	umor	Tumor	Total
ر	Non-Tumor	318	8	70	388
	Tumor	15		310	325
	Total	333		380	713
	Sensitivity: 0.934 Specificity: 0.815 Precision: 0.819 Accuracy: 0.881				

/lethod >	WBC Sequencing	Model Prediction
0.6% VAF	92.8% (64/69)	100% (69/69)
0.6% VAF	42.5% (74/174)	100% (174/174)



uniformity in variant prevalence across solid tumor cancer types in the Guardant plasma database.

Our bioinformatic model exhibits high sensitivity and specificity with WBC for discriminating tumor and non-tumor using only cfDNA. Our bioinformatic model has improved sensitivity for identifying non-tumor variants over WBC sequencing at low VAFs (<0.6%). In a paired plasma and WBC late stage cancer cohort, the majority of non-tumor variants were in known clonal hematopoiesis genes and variants of uncertain significance. No clinically actionable variants, except in ATM and CHEK2, were confirmed or annotated as non-tumor.

Conclusions

GUARDANTINFINGY

Cancer Res (2018) Jaiswal et al. Age-related clonal hematopoiesis associated with adverse outcomes. N Engl J Med (2014).